


Balance Software, EC 03



Presentation Screen

Print out

Balance Software EC 03

This version of the "Balance" software, can be connected to Sartorius and Mettler balances.

Calculations of Density, Weight change and Volume change.

For weight change and volume change the test "end time" is calculated from the start date and duration of the test. A reminder note can also be printed.

A material profile is created for each material with limits and target value. The data is stored per material and new measurements are added to existing material profiles.

The software collects the values from the connected instrument and plots a trend graph together with statistics. The results can be printed per material.

